

## TP8L1 rabbit pAb antibody

| Catalog No : | Source: | Concentration : | Mol.Wt. (Da): |
|--------------|---------|-----------------|---------------|
| A22768       | Rabbit  | 1 mg/ml         |               |

|                              |   |
|------------------------------|---|
| <b>Applications</b>          | WB,ELISA  |
| <b>Reactivity</b>            | Human,Mouse   |
| <b>Dilution</b>              | WB 1:500-2000 ELISA 1:5000-20000  |
| <b>Storage</b>               | -20°C/1 year  |
| <b>Specificity</b>           | TP8L1 Polyclonal Antibody detects endogenous levels of protein.   |
| <b>Source / Purification</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                         |
| <b>Immunogen</b>             | Synthesized peptide derived from human protein . at AA range: 40-120  |
| <b>Uniprot No</b>            | Q8WVP5  |
| <b>Alternative names</b>     |   |
| <b>Form</b>                  | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.  |
| <b>Clonality</b>             | Polyclonal  |
| <b>Isotype</b>               | IgG   |
| <b>Conjugation</b>           |   |
| <b>Background</b>            | similarity:Belongs to the TNFAIP8 family.,  |
| <b>Other</b>                 | TNFAIP8L1, Tumor necrosis factor alpha-induced protein 8-like protein 1 (TNF alpha-induced protein 8-like protein 1) (TNFAIP8-like protein 1) |

### Product Images:

#### Application Key:

WB-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation  
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

#### Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster

---

X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All

Species Expected

**Trademarks**

*All product names and trademarks are the property of their respective owners.*

**Regulatory Disclaimer**

*For life science research only. Not for use in diagnostic procedures.*

---

**Contact and Support:**

*To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).*

*To call, write, fax, or email us, please visit [www.aabsci.cn](http://www.aabsci.cn), contact information will be displayed.*

---