

## Nano-Tag9 mouse mAb(4C1) antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (Da):
A42040	Mouse	1 mg/ml	

<b>Applications</b>	WB
<b>Reactivity</b>	Species independent
<b>Dilution</b>	WB 1:2000-5000
<b>Storage</b>	-20°C/1 year
<b>Specificity</b>	The antibody detects the protein which have Nano-Tag9
<b>Source / Purification</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen</b>	Synthetic Peptide of Nano-Tag9
<b>Uniprot No</b>	
<b>Alternative names</b>	
<b>Form</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	
<b>Conjugation</b>	
<b>Background</b>	The Nano-tag is a new streptavidin-binding peptide for both the purification and the detection of Nano-tagged proteins. This peptide possesses nanomolar-affinity for streptavidin and therefore is termed Nano-tag. The nano-tags have two types, Nano-tag15(MDVEAWLGARVPLVET) and Nano-tag9(MDVEAWLGAR), which bind to streptavidin with dissociation constants of 4 nM and 17 nM, respectively.

### Other

### Product Images:

#### Application Key:

WB-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation

IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster

X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All

Species Expected

**Trademarks**

*All product names and trademarks are the property of their respective owners.*

**Regulatory Disclaimer**

*For life science research only. Not for use in diagnostic procedures.*

**Contact and Support:**

*To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).*

*To call, write, fax, or email us, please visit [www.aabsci.cn](http://www.aabsci.cn), contact information will be displayed.*